



# VETSCAN® SA & UA SAMPLE HANDLING

## BEST PRACTICE GUIDE

### URINE SAMPLE COLLECTION METHODS

#### Cystocentesis




Insert a small sterile syringe with needle into the bladder through the abdominal wall to aspirate urine.



Ideal method for collecting samples, requiring at least 2 veterinary professionals.

#### Catheter



Place a urinary catheter within the urethra and advance it into the bladder using a sterile syringe to collect the urine.



Requires at least 2 veterinary professionals and often, sedation of the animal.

#### Free Catch



Use a sterile collection container and collect a midstream sample.



Most at-risk for contamination and can be difficult to obtain.

### SAMPLE SIZE

**A benefit of in-clinic automated urinalysis is the small sample size volume required.**

- Acquire at least 1 mL of urine sample
- The SA requires approximately 650 µL (0.65 mL) of sample
- If running a UA test, additional sample may be required for the UA test strip

### SAMPLE PREPARATION

The SA does not centrifuge samples or use centrifuged samples as this process can destroy cellularity, casts, or crystals.



Invert tube to mix thoroughly 5 times immediately prior to running the sample.



Samples may be refrigerated in sealed, sterile containers for up to 4 hours.<sup>3</sup> Cold samples should be warmed to room temperature before testing. Refrigeration can enhance crystal formation in the urine.<sup>4</sup>



Urine samples may be tested immediately after collection for up to 1 hour post-collection, if stored at room temperature.<sup>1,2</sup>



Fresh, room temperature, well-mixed samples are ideal.

### DILUTION

#### When to dilute?

Dilution may be necessary when sediment elements overlap and the analyzer cannot identify individual elements.


#### What should be reviewed before dilution?

Review urine color, turbidity, analyzer warnings such as 'TNTC - dilute and retest' message, and 96 images to determine appropriate next steps.

#### How to perform the dilution?

If the provided images do not deliver clinically actionable information, a 1:2 dilution should be performed initially with increasing dilutions (1:4, 1:8) performed until images and/or quantitative or semi-quantitative values provide clinical direction.<sup>5</sup>


Dilution Ratio	Urine Volume	Sterile Saline Volume (0.9% NaCl)
1:2	0.4 mL	0.4 mL
1:4	0.2 mL	0.6 mL
1:8	0.1 mL	0.7 mL



Diluting a sample will alter its general properties. Only dilute the urine sample if necessary.

When performing dilution, it is important to:

- Use sterile 0.9% NaCl to dilute urine samples
- Invert tube to mix thoroughly 5 times immediately before running the sample



Contact Technical Support for additional information: ☎ 800 822 2947 @ DxSupport@zoetis.com

**Reference:** 1. Chew, Dennis and DiBartola, Stephen, Interpretation of Canine and Feline Urinalysis. Nestle Purina, Wilmington, De, 2004 2. Sink CA and Weinstein NM. Specimen Procurement In: Practical Veterinary Urinalysis. Ames, IA: John Wiley & Sons Inc. 2012. pgs. 12-15. 3. Sink CA and Feldman BF. Specimen Collection and Dipstick Analysis In: Laboratory Urinalysis and Hematology for the Small Animal Practitioner. Jackson, WY: Teton NewMedia. 2004. pgs. 3-18. 4. Sink CA, Weinstein NM. Routine Urinalysis: Microscopic Elements. Practical Veterinary Urinalysis. Chapter 5. Ames, IA: John Wiley & Sons Inc. 2012. 55-112. 5. Zoetis Data on File. Study No. TI-04858  
All trademarks are property of Zoetis Services LLC or a related company or a licensor unless otherwise noted. ©2020 Zoetis Services LLC. All rights reserved. VTS-00110